



TCA Threshold Testing in Red Wine

Determine your sensory threshold of TCA (2,4,6-Trichloroanisole) at Enartis Vinqury's neutral-environment sensory laboratory. During each 45 minute testing session, participants will be individually evaluated for their lowest level sensory threshold of TCA in a neutral red wine. Each person will perform a series of triangle tests using red wine containing 1-10 parts per trillion TCA. If the participant's threshold is determined to be in this range, they will receive a certificate specifying their sensory threshold level. Quantifying TCA levels is important to identify individuals that are able to recognize the defect. Winemakers, cellar masters and production staff can be used for quality control in a cork screening program or to identify TCA in the winery. Tasting room staff can also benefit by identifying corked bottles in the tasting room before they are served to consumers.

Testing must be done in an odor free environment, therefore please do not wear perfume or any other strong odors. The test will take approximately 30 minutes.

Reserve your space now!

There are a limited number of test times available. Please call (707) 838-6312 or fax this form back to (707) 838-1765 to set-up a test time.

Location

Enartis Vinqury Napa
2025 Redwood Road
Suite 9
Napa, CA 94558

Test times will be scheduled throughout the day. When registering please list your preferred test time:

December 16th _____

Fee: \$75 per person - *No refunds will be given the day of the seminar or thereafter.*

Name(s) _____

Company _____

Phone _____ Email/Fax* _____

Payment Method:

- Enclosed is a check payable to Vinqury
- Vinqury Account
- Visa/MasterCard/American Express

Account #: _____ Exp: _____ V-Code: _____

Cardholder: _____ Zip Code: _____

**If you do not receive confirmation via fax or email within 48 hours of submitting this form, please call to ensure that your space has been reserved.*